

TS-1064	Replacing the Peel Plate Tension Roller	Document Owner: Service Department	Effective Date: AUG-23-2022	Page 1 of 4
----------------	--	--	---------------------------------------	--------------------

Purpose

This guide provides the steps required to replace the peel plate tension roller.

Intended Use

If the peel plate tension roller needs to be replaced.

Required Material and Equipment

- 2.5mm allen wrench
- 1.5mm allen wrench

Intended Audience

Superuser

Prerequisites

N/A

Procedure

1. Remove the peel plate from the instrument. Loosen the single screw using a 2.5mm allen wrench, and disconnect the peel blow air line.



Figure 1.

2. Loosen the peel plate tension roller set screw found on the backside of the peel plate using a 1.5mm allen wrench.

TS-1064	Replacing the Peel Plate Tension Roller	Document Owner: Service Department	Effective Date: AUG-23-2022	Page 2 of 4
----------------	--	--	---------------------------------------	--------------------



Figure 2.

3. Remove the faulty shaft, by pulling up.



Figure 3.

4. Place the new shaft into the open position with the push nut on top as shown in Figure 4.

TS-1064	Replacing the Peel Plate Tension Roller	Document Owner: Service Department	Effective Date: AUG-23-2022	Page 3 of 4
----------------	--	--	---------------------------------------	--------------------



Figure 4.

5. Tighten the set screw.



Figure 5.

6. Install the peel plate on the instrument.

TS-1064	Replacing the Peel Plate Tension Roller	Document Owner: Service Department	Effective Date: AUG-23-2022	Page 4 of 4
----------------	--	--	---------------------------------------	--------------------



Figure 6.

If you need additional assistance, please call our service department at 314-298-9800 or email service@scinomix.com.

Definitions

N/A

Revision History

Version:	Change:	Effective Date:	Approved by:
A	NEW	MAY-04-2021	B. Petti
B	Content review Added <i>Intended Use</i>	AUG-23-2022	Service