

<b>TS-1064</b>	<b>Replacing the Peel Plate Tension Roller</b>	<b>Document Owner:</b> Service Department	<b>Effective Date:</b> May-04-2021	<b>Page 1 of 4</b>
----------------	--	--	---------------------------------------	--------------------

## Purpose

This guide provides the steps required to replace the peel plate tension roller.

## Required Material and Equipment

- 2.5mm allen wrench
- 1.5mm allen wrench

## Intended Audience

Field Service Representative from Scinomix or approved vendors.

## Prerequisites

N/A

## Procedure

1. Remove the peel plate from the instrument. Loosen the single screw using a 2.5mm allen wrench, and disconnect the peel blow air line.



Figure 1.

2. Loosen the peel plate tension roller set screw found on the backside of the peel plate using a 1.5mm allen wrench.

<b>TS-1064</b>	<b>Replacing the Peel Plate Tension Roller</b>	<b>Document Owner:</b> Service Department	<b>Effective Date:</b> May-04-2021	<b>Page 2 of 4</b>
----------------	--	--	---------------------------------------	--------------------



Figure 2.

3. Remove the faulty shaft, by pulling up.



Figure 3.

<b>TS-1064</b>	<b>Replacing the Peel Plate Tension Roller</b>	<b>Document Owner:</b> Service Department	<b>Effective Date:</b> May-04-2021	<b>Page 3 of 4</b>
----------------	--	--	---------------------------------------	--------------------

- Place the new shaft into the open position with the push nut on top as shown in Figure 4.



Figure 4.

- Tighten the set screw.



Figure 5.

<b>TS-1064</b>	<b>Replacing the Peel Plate Tension Roller</b>	<b>Document Owner:</b> Service Department	<b>Effective Date:</b> May-04-2021	<b>Page 4 of 4</b>
----------------	--	--	---------------------------------------	--------------------

6. Install the peel plate on the instrument.



Figure 6.

*If you need additional assistance, please call our service department at 314-298-9800 or email [service@scinomix.com](mailto:service@scinomix.com).*

## Definitions

N/A

## Revision History

<b>Version:</b>	<b>Change:</b>	<b>Effective Date:</b>	<b>Approved by:</b>
A	NEW	MAY-04-2021	B. Petti